

Subject: Atomic Absorption Spectrophotometer – PR4922785

REQUEST FOR PRICE QUOTATION

Dear Sir/Madam,

We kindly request you to submit a price quotation for

1 ea. Bench Top Atomic Absorption Spectrophotometer – Specifications attached.

The quotation may be for **Atomic Absorption Spectrophotometer (Graphite Furnace)** or for an alternate brand/make/model that includes the following **per “attached” specification listing** salient features/specifications:

Attachment: Specification of Atomic Absorption Spectrophotometer (Graphite Furnace)

Your quotation shall include FOB delivery to **State site: New Jersey, NY, in Europe: Antwerp/Belgium, final destination: Yerevan, Armenia** and the following:

- (a) Documentation/confirmation for product that is being offered, that the item meets the essential salient characteristics according to the attached specification.
- (b) Price(s) / cost of item(s) including packing and delivery charges.
- (c) Delivery terms and conditions.
- (d) Validity period of the quotation.
- (e) Any other information which may be required for the particular item quoted, e.g., warranty

Please submit your quotation no later than December 28, 2015 – 16:00 hours (German local time) to:

American Consulate General Frankfurt, Germany
Regional Procurement Support Office (RPSO)
Procurement Agent: Hannelore Godecki
Giessener Strasse 30
60435 Frankfurt am Main – Germany

or email to: GodeckiH@state.gov

Sincerely,

[Hannelore Godecki]

Dated 12/17/2015 HG.

**Specification for Procurement of
Atomic Absorption Spectrophotometer (Graphite Furnace)**

Model	<ul style="list-style-type: none"> • Bench Top
Lamp mount	<ul style="list-style-type: none"> • Minimum six or more lamps with facility to accommodate high efficiency/energy lamps with built in power supplies. • Lamp selection, alignments and operating current should be software controlled.
System	<ul style="list-style-type: none"> • Incorporate preferably a true double beam flame spectrometer • The sample and reference beams are measured simultaneously for enhanced precision and detection limits • Zeeman Graphite furnace atomizer in an integrated manner
Monochromator	<ul style="list-style-type: none"> • For automatic wavelength selection • Holographic grating with minimum of 1800 lines/mm • Wavelength range 195-850 nm
Automatic slit width	<ul style="list-style-type: none"> • 0.2, 0.7 and 2.0 nm is preferred
Detectors	<ul style="list-style-type: none"> • Photomultiplier • Wavelength length 195-850 nm
Background Correction	<ul style="list-style-type: none"> • High speed deuterium background correction for flame operation and Longitudinal Zeeman-effect background correction • The magnet is automatically switched on during the atomization step
Lamps	<ul style="list-style-type: none"> • Wireless Coded Hollow Cathode Single element lamps for the elements viz. Na, K, Ca, Mg, Cr, Mn, Fe, Ni, Cu, Zn, Cd, Al, Si, P, Co, Sb, Ba, Be, Mo, Se, Sr • High Efficiency lamps of single element for Pb, Hg, As. • Suitable standards for all the elements need to be quote. • Detection limits with best possible instrument configurations should specify for the above elements
Flame Atomizer	<ul style="list-style-type: none"> • Software controlled oxidant and fuel monitoring • Remote ignition of flame and safe switching between fuel gases • Gas control must be programmable through software for every element in flame mode. • Safety interlocking for burner, drain and gases • Software controlled automatic alignment of flame in the light beam • Inert spray chamber and corrosion resistant nebulizer • Single slot titanium burner heads of 10 length suitable for air/acetylene and 5cm length for nitrous oxide
Graphite Furnace Atomizer	<ul style="list-style-type: none"> • Fully computer controlled • Uniform heating of Graphite Tube for better and consistent results • Transversely Heated Graphite Atomizer with built in integrated platform for all the elements • Temperature program: ambient to 2300°C in a step of 10°C or better, minimum 12 program

	<ul style="list-style-type: none"> • Variable gas flow for furnace for low minimum consumption • an external protective gas stream around the graphite tube prevents the entrance of outside air to maximize tube life • An internal purge gas goes through the graphite tube to remove the volatilized matrix vapors during and thermal pretreatment • The two gas streams are computer controlled independently • Pneumatically operated software controlled furnace opening and closing • Argon as inert gas with minimum gas consumption • Imported Self-priming water recirculation cooling systems • Furnace auto sampler with minimum 50 samples, references and pipette wars containers • Control of auto sampler through spectrometer software • Sample and Reagent with Wide Volume, Selectable in increments of 1 µl
Furnace Viewing Color Camera	<ul style="list-style-type: none"> • For taking images of the events taking place inside the graphite furnace cuvette for easy sample viewing, optimization and method development.
Software/System	<ul style="list-style-type: none"> • Licensed control, application, data management and complete quality control software package, compatible with windows 7 operating system with their respective license keys. • Branded desktop computer preloaded with latest operating system and AAS control/application software and laser printer.
Accessories	<ul style="list-style-type: none"> • Suitable Imported air compressor for AAS • Fully automated continuous Flow Injection system for mercury and hydride-forming elements. (Should be coupled with Flame and Furnace technique to achieve the lowest possible detection limit) Elements Detection limit like As - 0.005ppb, Se – 0.005ppb and Bi – 0.006ppb. • Two numbers of Cylinders (with safety certificate) for all the required gases • with suitable two stage regulators, wall mountable filters and brackets and tubing for installation • Required exhaust/fume hood system for AAs – Supply and Installation. • Required tubing, connecting cables and fittings for installation and commissioning of AAS and accessories • Supply of Graphite tubes – quantity 5
Warranty	<ul style="list-style-type: none"> • Minimum 2 Years